

Workforce Safety & Insurance
 1600 E Century Ave Ste 1
 PO Box 5585
 Bismarck, ND 58506-5585
 701.328.3800
 1.800.777.5033
 www.WorkforceSafety.com

Revised Document Date: 07/20/2015

**Pharmacy Benefit Management (PBM) Program
 FORMULARY/PRODUCT RESTRICTIONS**

Key:

B Available as brand name only
G Available as a generic equivalent

\$ Average cost < \$50 per prescription
\$\$ Average cost \$50 to \$99 per prescription
\$\$\$ Average cost \$100 to \$249 per prescription
\$\$\$\$ Average cost \$250 to \$499 per prescription
\$\$\$\$\$ Average cost > \$499 per prescription

F Formulary: Indicates that the product is available without prior authorization.
N/F Non-Formulary: Indicates that the product is not available. An alternative F (Formulary) or PA (Prior Authorization) agent must be selected.
PA Prior Authorization: Indicates that the product require prior approval by WSI before being dispensed to the injured worker.

| 51 Digestive Aids | Description | Brand Name | B/G | Cost | Status |
|---------------------------------------|--|---------------------|-----|------|--------|
| <i>5120 Digestive Enzymes</i> | | | | | |
| | Lactase Enzyme Cap 250mg | Lactrase | B | | PA |
| | Lactase Enzyme Tab 3000 units | Lactaid | G | | PA |
| | Lactase Enzyme Tab 9000 units | Lactaid Fast Acting | G | | PA |
| | Lactase Enzyme Tab Chew 9000 units | Lactaid Fast Act | G | | PA |
| | Sacrosidase Soln 8500 IU/ml | Sucraid | B | | PA |
| <i>5199 Digestive Aids - Mixtures</i> | | | | | |
| | Amylase-Lipase-Protease Cap Del Rel 15,000u/3000u/9500u | Creon | B | | PA |

| | | | |
|---|----------------------|---|----|
| Amylase-Lipase-Protease Cap Del Rel 16,000u/3000u/10,000u | Zenpep | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 17,500u/4200u/10,000u | Pancreaze | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 20,000u/4500u/25,000 | Ultrase | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 25,000u/4000u/25,000u | Pancrecarb MS-4 | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 27,000u/5000u/17,000u | Pancrelipase, Zenpep | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 27,6000u/13,800u/27,600u | Ultresa | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 30,000u/6000u/19,000u | Creon | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 30,250u/8000u/28,750u | Pertzye | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 39,000u/12,000u/39,000u | Ultrase MT12 | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 40,000u/8000u/45,000u | Pancrecarb MS-8 | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 41,400u/20,700u/41,400u | Ultresa | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 43,750u/10,500u/25,000u | Pancreaze | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 46,000u/23,000u/46,000u | Ultresa | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 52,000u/16,000u/52,000u | Pancrecarb MS-16 | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 55,000u/10,000u/34,000u | Zenpep | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 58,500u/18,000u/58,500u | Ultrase MT18 | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 60,000u/12,000u/38,000u | Creon | B | PA |
| Amylase-Lipase-Protease Cap Del Rel 65,000u/16,000u/57,500u | Pertzye | B | PA |

| | | | | |
|--|-----------------|---|------------|----|
| Amylase-Lipase-Protease Cap Del Rel 61,000u/21,000u/37,000u | Pancreaze | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 65,000u/20,000u/65,000u | Ultracaps MT-20 | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 70,000u/16,800u/40,000u | Pancreaze | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 82,000u/15,000u/51,000u | Zenpep | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 109,000u/20,000u/68,000u | Zenpep | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 120,000u/24,000u/76,000u | Creon | B | \$\$\$\$\$ | PA |
| Amylase-Lipase-Protease Cap Del Rel 136,000/25,000/85,000u | Zenpep | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 180,000u/36,000u/114,000u | Creon | B | | PA |
| Amylase-Lipase-Protease Cap Del Rel 218,000u/40,000u/136,000u | Zenpep | B | | PA |
| Amylase-Lipase-Protease Tab 39,500u/10,400u/39,150u | Viokase 10 | B | | PA |
| Amylase-Lipase-Protease Tab 78,300u/20,880u/78,300u | Viokase 20 | B | | PA |